

Překlad z anglického jazyka

XIANTAO ZHONGTAI PROTECTIVE PRODUCTS CO., LTD

Osvědčení o analýze

Název výrobku	Zesílený chirurgický plášť		
Model	L (134x147 cm)	Datum odebrání vzorku	10.1. 2021
Datum výroby	5.1. 2021	Platí do	01. 2024
Č. šarže	ZT-201225-1	Datum zprávy	15.1. 2021
Množství	100 ks	Norma	EN 13795-1:2019
Obsah analýzy	Požadavek		Výsledek analýzy
Požadavek týkající se materiálů	Hlavní materiál	SMS materiál, 40 g/m ² , modrý PE 40 g/m ² , bílý, s PP potažením	Vyhovuje
	Tvar	Otvor vzadu, úpletové manžety, se suchým zipem u krku, v pase 4 úvazy. Zesílení v přední části a na rukávech	Vyhovuje
	Velikost	Délka 134+/- 2 cm, šířka 147+/-2cm	Vyhovuje
	Technologie	Zatavení ultrazvukem	
	Značení (štítek)	Vytisknuto na sáčku z kombinovaného materiálu papír - plast	Vyhovuje
	Balení	1 ks plášť a 2 ks ručníků v sáčku z kombinovaného materiálu papír – plast, netkaná textilie, baleno 30 ks v kartonu	Vyhovuje
	Karton	51x42,5x32,5 cm, pětivrstvý karton, potištěný	Vyhovuje

Závěr: Výše uvedené výrobky splňují normu EN 13795-1:2019.

Kontroloval:

Analýzu provedl:

Firemní razítko:

XIANTAO ZHONGTAI PROTECTIVE PRODUCTS CO., LTD

XIANTAO ZHONGTAI PROTECTIVE PRODUCTS CO., LTD.

Certificate of Analysis

Name of Production	Reinforced Surgical Gown		
Model	L (134x147 cm)	Sampling Date	2021.01.10
Manufacture Date	2021.01.05	Expiry Date	2024.01
Lot No	ZT-201225-1	Data of Report	2021.01.15
Quantity	100pcs	Standard	EN 13795-1: 2019
Analysis Content	Requirement		Analysis Result
Physical Requirement	Main Material	The material SMS is 40gsm, blue PP coated PE 40gsm, white	Confirms
	Shape	Back opening, knitted cuffs, with velcro at neck, four ties at waist Reinforced at front chest and sleeves	Confirms
	Size	Length 134+-2cm, width 147+-2cm	Confirms
	Technology	Ultrasonic Sealing	Confirms
	Label	Printing on paper-plastic bag	Confirms
	Packing	1pc Gown and 2pcs tissues in a paper-plastic bag, non-woven wrapped 30pcs in a carton	Confirms
	Carton	51x42.5x32.5cm, five-layer carton, printed	Confirms

Conclusion: The above products meet the standard EN 13795-1: 2019

Checker:

廖长群

Analyst:

胡雨舟



扫描全能王 创建