

### Osvědčení o analýze

Název výrobku	Třívrstvá ústenka s gumičkou		
Model	17,5 × 9,5 cm	Datum vzorkování	25. 10. 2023
Datum výroby	20. 10. 2023	Datum expirace	10. 2028
Číslo šarže	231020	Datum zprávy	27. 10. 2023
Množství	200 ks	Norma	EN 14683:2019+AC:2019 Typ IIR
Obsah analýzy	Požadavek		Výsledek analýzy
Fyzické požadavky	Hlavní materiál	Materiál tvoří 2 vrstvy netkané textilie z polypropylenu + 1 vrstva netkaná textilie z MB, netoxické, bez zápachu, hypoalergenní, prodyšné, hydrofobní, bez skelných vláken	Vyhovuje
	Tvar	Obdélník se třemi sklady, tvarovatelná fixační část v oblasti nosu a úst	Vyhovuje
	Účinnost filtrace bakterií	Větší nebo rovno 95 %	Vyhovuje
	Poutka za uši	Gumičky ze spandexového vlákna	Vyhovuje
	Technologie	Ultrazvukové tepelné utěsnění	Vyhovuje
	Rozměry	Délka 17,5+0,2 cm; šířka 9,5+0,2 cm	Vyhovuje
	Balení	50 ks v jedné krabičce, 2000 ks v kartonu	Vyhovuje
	Krabička	18 × 10 × 7 cm, s potiskem	Vyhovuje
	Karton	52 × 38 × 30 cm, pětivrstvý karton s potiskem	Vyhovuje

Závěr: Výše uvedené výrobky splňují normu EN 14683:2019+AC:2019 Typ IIR

Zkoušející: Xu Chen

Analytik: Zhang Keping

### Certificate of Analysis

Name of Production	3 layer Face Mask with earloop		
Model	17.5*9.5cm	Sampling Date	2023.10.25
Manufacture Date	2023.10.20	Expiry Date	2028.10
Lot No	231020	Data of Report	2023.10.27
Quantity	200pcs	Standard	EN 14683:2019+AC:2019 Type IIR
Analysis Content	Requirement		Analysis Result
Physical Requirement	Main Material	The material is 2ply pp nonwoven+1ply MB nonwoven, non-toxic, odorless, Hypoallergenic, breathable, hydrophobic, free of glass fibers	Conform
	Shape	Rectangle with three folding, a formable fixation part in the area of the nose and mouth	Conform
	Bacterial filtration efficiency	Greater than or equal to 95%	Conform
	Ear-loop	Elastic Band made of Spandex filament	Conform
	Technology	Ultrasonic heat-seal	Conform
	Size	Length 17.5+/-0.2cm; Width 9.5+/-0.2cm	Conform
	Packing	50pcs in one box, 2000pcs in a carton	Conform
	Box	18x10x7cm, printed	Conform
Carton	52x38x30cm, five-layer carton, printed	Conform	

Conclusion: The above products meet the standard EN 14683:2019+AC:2019 Type IIR

Checker: Xu Chen

Analyst: Zhang Keping